



Ciljna skupina

Prezbiopi koji većinu dana provode u zatvorenom prostoru, radeći za laptopom, čitajući razna dokumenta ili u razgovorima sklijentima.
Životna dob 45+.

Bilješka: Strogo je zabranjena uporaba ove vrste leće u prometu.

Tehnologija

OPTOTECH OFFICE 16

Način izrade

Unutarnja izrada digresivne površine.

Poručivanje dioptrija

1/4

Progressivna duljina

18 mm

Minimalna visina ugradnje

20 mm

Način poručivanja

Dioptrija za blizinu

Odabir digresije na osnovu tablice (dolje) ili po želji

Koristiti razmak zjenica za daljinu, po mogućnosti za svako oko pojedinačno (IPD).

Optimizacija

Specijalno tanjenje za leću u plus rasponu na osnovu oblika okvira i parametra za fitovanje leće u okvir. Ako se leća poručuje bez opcije Optimizacija, topografija leće bit će decentrirana nazalno u odnosu na geometrijsko središte za 2,5 mm.

Ocjene korisnika

Vid – daljina



Vid – blizina



Spontano privikavanje



Vid – srednja daljina



Udobnost tijekom korištenja



Korištenje elektroničkih uređaja




Oznake na lećama



SELEKCIJA DIGRESIVNE JAČINE

Adicija	Udaljenost do 1 m	Udaljenost do 2 m	Udaljenost do 3 m
1,00	-	-	0,75
1,25	-	0,75	1,00
1,50	0,75	1,00	1,25
1,75	1,00	1,25	1,50
2,00	1,25	1,50	1,75
2,25	1,50	1,75	2,00
2,50	1,75	2,00	2,25
2,75	2,00	2,25	-
3,00	2,25	-	-



MATERIJALI

1.50 CR-39	1
1.53 Trivex	2
1.56 Mid. index	3
1.59 Poly	4
1.60 High Index	5
1.60 Tribrid	6
1.67 High Index	7
1.74 High Index	8

DORADE

Bez dorade	W	Transitions GEN 8 Sapphire	F	Nupolar Gray	S
Remuve 400	H	Transitions GEN 8 Amethyst	M	Nupolar Gray Light	L
Remuve 420	R	Transitions GEN 8 Amber	A	Nupolar Green	E
Remuve 420 Photomatic Brown	O	Transitions GEN 8 Emerald	E	Transitions XTRActive Polarised Gray	X
Remuve 420 Photomatic Gray	U	Transitions XTRActive Brown	C	Photomatic Brown	B
Transitions GEN 8 Brown	Z	Transitions XTRActive Gray	L	Photomatic Gray	I
Transitions GEN 8 Gray	N	Transitions XTRActive Green	K	DriveWear	D
Transitions GEN 8 Green	V	Nupolar Brown	B		